

# Lab Blender Model MD-326S





No.02016326 (Zhejiang, CHINA)

# ...Quality is Everything...

**Telephone:** +86-579-87537698, +86-13819904200 **Fax**: +86-579-8753 696 Rev.B

Website: <a href="http://www.rigchina.com">http://www.rigchina.com</a> Email: <a href="mailto:sales@rigchina.com">sales@rigchina.com</a> MSN:rigchina@hotmail.com</a> SKYPE:

Address: No.80-82, Qiude Road, West Cheng Industrial Estate, Yongkang city Zhejiang Province, China SKYPE:rigchina



## **LABORATORY BLENDER MD-326S**

## 10 levels of variable speed

The Model MD-326S Lab Blender with 10 levels of variable speed, making easier for adjusting and choosing speed level.

3800.7600.11400.15200.19000.22800.26600.30400.34200.38000

## 3 auto-timer program

The Model MD-326S Lab Blender provides 3 auto-timer: 30, 60, and 120 seconds.





## Durable

BPA Free blender jar, made of #304 stainless steel, available for cold and hot temperature. Full metal upper and lower clutch.







## **Quiet Sound**

The Model MD-326S Lab Blender brings a new invention of the new concept of air-flowing system to reduce the heat of the motor. These features allow it to operate in lower noise.

High Speed Lab Blender Specifications			
Part No.	Model	Description &Technical Parameters	
326S	MD-326S	10 Preset of Speeds:	
		3800,7600,11400,15200,19000,22800,26600,30400,34200,38000	
		RPM	
		Variable Speeds Range: 0-38,000 RPM.	
		220V±5%AC; 50/60Hz; Power: 1800W;	
		Capacity:1800ml;	
		3 auto-timer: 30, 60, and 120 seconds.	
		Net weight /Gross Weight: 6Kg / 10 kg	



# **Warranty and Returns**

## Warranty

Rigchina Group Company warrants its products to be free from defects in material and workmanship for a period of 12 months from the time of shipment. If repair or adjustment is necessary, and has not been the result of abuse or misuse within the twelve-month period, please return, freight prepaid, and correction of the defect will be made without charge.

Out of warranty products will be repaired for a nominal charge.

Please refer to the accompanying warranty statement enclosed with the product.

## Returns

For your protection, items being returned must be carefully packed to prevent damage in shipment and insured against possible damage or loss. Rigchina will not be responsible for damage resulting from careless or insufficient packing.

Before returning items for any reason, authorization must be obtained from Rigchina Instrument Company. When applying for authorization, please include information regarding the reason the items are to be returned.

#### Our shipping address is:

#### RIGCHINA GROUP COMPANY

Full Address: No.80-82, Qiude Road, West Cheng Industrial Estate,

Yongkang city, Zhejiang Province,

China

Postcode: 321300

Telephone: 0086-579-87537698

Fax: 0086-579-8753 696 Email: service@rigchina.com

Website: <a href="http://www.rigchina.com">http://www.rigchina.com</a>